

WEBINAR

Optimising Spirometry Endpoints in Clinical Trials



Grant Sowman

VP Clinical Services
Vitalograph



Holly Wilson

Pulmonary Function Technologist
RPFT

Why this matters

- Spirometry drives endpoints
e.g. FEV1: FVC: % Predicted
- Missing data and variability reduces statistical power
- ATS-ERS 2019 Spirometry testing guidelines - does the approach need to be disease specific?



The hidden risk

- Are we automatically rejecting good FEV1 and FVC data that does not meet a specific ATS-ERS QA criteria due to a patient's physiology e.g. IPF (restrictive) vs COPD (obstructive)?

- Are we creating bias from these exclusions?

- Are we unnecessarily doing repeat visits?



Learning objectives

1 Interpret ATS-ERS 2019 Test QA criteria

2 Evaluate disease-specific risk

3 Assess limits of uniform QA

4 Apply disease-specific QA



ATS-ERS objective QA test criteria errors explained



Excessive Back Extrapolation

- BEV in excess of 5% of FVC or 0.10L



Submaximal Inhalation

- Failing to reach TLC prior to forced manoeuvre.
- Submaximal effort.
- Lack of test understanding.



Cough in 1st Second

- An audible cough that raises suspicion of incorrect expiratory technique.
- A cough that is visible on the graphs which impacts the forced expiratory flow in any way.



EOFE Criteria Not Met

- Early cessation of expiration
- Failing to achieve a 1s plateau
- Failing to achieve an FVC within the repeatability criteria of other FVCs within the session.



Glottis Closure

- Visible and abrupt fluttering of expiratory flow and volume on graphs which ends the manoeuvre or is followed by a second expiratory trace.

Example impact on endpoints

ISSUE	IMPACT ON FEV1	IMPACT ON FVC
Slow Start (BEV)	↓ Reduced	↓ Reduced
Glottis Closure	–	↓ Reduced
Cough (1st second)	↑ Increased	–
Submaximal Inhalation	↓ Reduced	↓ Reduced
EOFE	–	↓ Reduced

ATS-ERS subjective QA test criteria errors explained




Poor Coordination

- Participant struggling to perform relaxed tidal breathing following by a forced manoeuvre.



Artefact

- Any significant deviation in flow caused by cough, mouthpiece obstruction, changes in effort.



Drift

- Flow-signal error due to flow-head exposure to external airflow while zeroing.



Submaximal Effort

- Lack of test understanding
- Low expiratory effort applied in part or throughout.



Leak



- Poor lipseal
- Equipment fault or poor assembly.

Example impact on endpoints

ISSUE	IMPACT ON FEV1	IMPACT ON FVC
Drift	↑ Artificially high	↑ Artificially high
Submaximal effort	↓ Reduced	↓ Reduced
Leak	↓ Reduced	↓ Reduced
Poor coordination	↓ Reduced	↓ Reduced
Artefact -	↑ or ↓ (High or Reduced)	↑ or ↓ (High or Reduced)



Test acceptability % by disease

Obstructive




Test acceptability % by disease

Neuromuscular and Restrictive




Comparison COPD vs IPF

Reasons for failed Spirometry tests **IPF**



Comparison COPD vs IPF

Reasons for failed Spirometry tests **COPD**



Example disease-specific factors in patients

Restrictive



e.g. IPF often struggle to meet the back extrapolated volume criteria due to lower FVC values.

Severe Obstructive



e.g. COPD there was a lot of artefact particularly at the end of the test.

Neuromuscular Disorder



Often struggle with start of test and keeping maximal effort throughout the manoeuvre.

Assessing limitations

- Are there any risks with an automated QA approach?
- Could restrictive physiology (e.g. IPF) impact a patient's ability to achieve acceptable data?

Definition of BEV (Slow-start)



Back Extrapolated Volume measures air exhaled before maximal effort in forced expiratory manoeuvres.

BEV Criterion Threshold




BEV must be $\leq 5\%$ of Forced Vital Capacity or 0.100 liters, whichever is greater

Acceptability and Usability



BEV requirement applies equally to acceptability and usability for both FEV₁ and FVC measurements.


What are the risks?



PARAMETER	TEST 1	TEST 2	TEST 3
FVC (L)	1.74	1.73	1.89
FEV1 (L)	1.37	1.37	1.44
PEF (L/min)	284	266	349
BEV (L)	0.11	0.08	0.11
TPEF (ms)	145	117	111
Repeatability - FVC (ml)	-140	-150	10
Repeatability - FEV1 (ml)	-70	-70	0
Start of Test Acceptable	NO	YES	NO

What are the risks?

Is this artefact/cough
consistent with the
automated review?





Clinical trial impact

- Could we be missing valid or usable data by automatically rejecting tests that don't meet a specific QA criteria without considering a patient unique physiology?
-
- Is there room to have a more flexible approach to apply an over-read of tests?



Applying a disease-specific approach to test QA

- Allow a human review of tests flagged as not meeting QA criteria that is disease specific.
-
- Flexible approach – document the scientific reasoning of why a test is acceptable despite a QA criteria not being met

Site-level changes

Tailored coaching for each patient population.

Adjust the time in-between tests to allow patient recovery.

Allow the Physiologist/ Investigator to over-ride any automated review with documented scientific reasoning.

Sponsor/CRO changes

Over-reading QA criteria

Allow for borderline test cases to be assessed by an expert.

Develop a practical spirometry QA framework

Aligned to your disease population.

Site training strategy

Adapt it to your patient population focus on the specific QA errors.





Human over-reader QA feedback

- Can identify any patient specific physiological patterns

- Determine if an automated QA flag is disease related or a true error.

- Provide clear instructions back to site staff conducting Spirometry testing how to improve coaching based on critical thinking and not a generic guide.

Key takeaways

- ✓ Understand the specific test QA challenges for the patient population in your clinical trial.
- ✓ Errors are disease-specific.
- ✓ Develop a disease-specific Spirometry QA criteria to be applied across all stakeholders, Sponsor, CRO, Investigator Site.
- ✓ Include an element of human over-reader review with a flexible approach to FEV1 and FVC test acceptability.



Thank you